





Introducing the game-changing elvatome® 2.0- the ultimate dental instrument that combines CUTTING-EDGE PRECISION and UNRIVALED STRENGTH to

revolutionize tooth extractions. With its spade-shaped sharp tip, elvatome® 2.0 effortlessly severs the periodontal ligament while expanding the alveolar bone, making extractions faster and easier than ever before.

Our NEW PATENT-PENDING TWIST

HANDLE promotes superior ergonomics, guiding clinicians to the recommended "Push and Twist" motion, ensuring precision and accuracy with every use. elvatome® 2.0 has a thicker construction than its predecessor, maximizing rigidity, and reducing clinician fatigue when working in interproximal areas.

Crafted with single rod construction and a weighted stainless steel stabilizer teardrop, elvatome[®] 2.0 maximizes strength and durability while effortlessly balancing the tip, allowing for easier entry into proximal areas with force.

This revolutionary instrument is a must-have in any dental practice's armamentarium, particularly for those seeking to Protect the Tissue, Preserve the Bone, and Enhance Patient Experience during tooth extractions.

Upgrade your dental practice with elvatome[®] 2.0 and experience a new level of **PRECISION**, **STRENGTH**, and **EASE** during every extraction.

Weighted Stabilizer Teardrop

TWIST Handle Patent-Pending

Single Rod Construction

Concave Sharp

Blade

38021.2 Straight Spade elvatome[®] 2.0

38026.2 3mm Curved elvatome[®] 2.0

Ideal for premolar to premolar in the upper or lower arch.

Ideal for premolar to premolar in the upper or lower arch.





38024.2 Left Angle Spade elvatome[®] 2.0 Ideal for mandibular molars and root tips.

38025.2 **Right Angle Spade** elvatome[®] 2.0 Ideal for mandibular molars and root tips.

Scan to See in Action



38022.2 **Distally Angled** Spade elvatome[®] 2.0

Ideal for distal, lingual, maxillary premolar and molars.



38023.2 **Mesially Angled** Spade elvatome® 2.0


Ideal for mesial, buccal-mesial, maxillary premolar and molars.



38100.2

elvatome[®] 2.0 Kit (5 Pack)

Kit includes 1 each of: Straight Spade, Distally Angled Spade, Mesially Angled Spade, Left Angled Spade, Right Angled Spade with Sterilizable Cassette.





For more information on our products and how to order visit **tbsdental.com.**

